

Empty EVA 2 Port Bags

Ethyl Vinyl Acetate (EVA) bags are constructed from a single-web, medical grade laminate of four different films, with a fluid contact surface of EVA, with two layers of low-density polyethylene (LDPE) enveloping a layer of ethyl vinyl alcohol (EVOH).

These bags are specifically designed to provide strong seals, extraordinary robustness, superior flex, crack and pin-hole resistance, plus desirable gas and moisture barrier performance.

- EVA single-use bags are designed for bio-pharmaceutical, laboratory and diagnostic applications
- 100% PVC, DEHP, and Latex-Free
- Practical and sturdy carrying handles
- Spike Port
- Injection Site
- Can be filled by gravity, with syringes or with single-station filling devices
- Transparent for easy admixture inspection
- Individual blister package rather than in case



WP-EVA250ML



WP-EVA100ML

Item	Description	Quantity
WP-EVA50ML	50 mL Empty EVA 2 Port Bags	50/Case
WP-EVA100ML	100 mL Empty EVA 2 Port Bags	50/Case
WP-EVA250ML	250 mL Empty EVA 2 Port Bags	50/Case

**CONTACT US TODAY TO EXPERIENCE
UNBEATABLE VALUE & SERVICE**

(800) 722-7865 | www.preferredmedical.com

Delivering Efficiency to Healthcare®

402 BNA Drive, Suite 500 / Nashville, TN 37217
www.preferredmedical.com / 615.722.7865