

Why Wolf-Pak® IV Start Kits?

These single-use Kits are designed to help medical staff to streamline their processes, thereby reducing waste and saving time. All the components of each tray are pre-packed in order of sequence, allowing for faster progression through each procedure. Our kits are Latex-Free and are 50/Case.



IV5500 includes:

- (1) Pair of Latex-Free Exam Gloves
- (1) Latex-Free Tourniquet
- (1) Alcohol Prep Pad
- (1) PVP Pad
- (2) Gauze Sponge, 2"x 2"
- (1) Transparent Dressing
- (1) Roll Tape
- (1) Dressing Change Label



IV5060LF includes:

- (1) Pair of Latex-Free Exam Gloves
- (1) Latex-Free Tourniquet
- (1) ChloraPrep™ Frepp™ 1.5 mL Applicator
- (1) Gauze Sponge, 2"x 2"
- (1) Opsite® Dressing
- (1) Roll Tape
- (1) Dressing Change Label



IV5090LF includes:

- (1) Pair of Latex-Free Exam Gloves
- (1) Latex-Free Tourniquet
- (1) Prevantics® Swab
- (2) Gauze Sponge, 2"x 2"
- (1) Transparent Dressing
- (1) Roll Tape
- (1) Dressing Change Label



IV6010LF includes:

- (1) Latex-Free Tourniquet
- (1) Hubguard Securment Device
- (1) Prevantics® Swab
- (2) Gauze Sponge, 2"x 2"
- (1) Transparent Dressing
- (1) Roll Tape
- (1) Dressing Change Label

**CONTACT US TODAY TO EXPERIENCE
UNBEATABLE VALUE & SERVICE**

(800) 722-7865 | www.preferredmedical.com

Delivering Efficiency to Healthcare®

402 BNA Drive, Suite 500 / Nashville, TN 37217
www.preferredmedical.com / 615.722.7865